

OptraSculpt® Pad

The modelling instrument for composite materials

Efficient
contouring
due to non-stick effect



OptraSculpt® Pad

The modelling instrument for composite materials

OptraSculpt® Pad is a contouring instrument with special foam pad attachments, which is designed for the efficient, non-stick forming and shaping of composites.



Shaping and contouring with OptraSculpt Pad



Shaping and contouring with a metal spatula

Non-stick shaping and contouring

The non-stick attachments of OptraSculpt Pad enable composite materials to be shaped and contoured with ease, without leaving any unwanted marks. Composite restorations with smooth and even surfaces are fabricated with utmost efficiency.




Result achieved with OptraSculpt Pad



Result achieved with a metal spatula

Creation of smooth and even surfaces

Due to the special material of the pads, natural-looking restorations are easily accomplished in only a few steps. The highly flexible synthetic foam pads optimally adjust to the anatomical contours and allow a smooth modelling.



Microscopic image of the foam




Reference scale 1



Reference scale 2


Professional esthetic results

The reference scales on the instrument handle assist in the creation of esthetic and anatomically-correct restorations. The markings allow the clinical situation to be compared with the ideal average tooth width proportions and angular alignments in the upper anterior dentition.



Benefits at a glance

- The foam pad attachments ensure a non-stick adaptation.
- The instrument does not leave any marks and therefore ensures the creation of even surfaces.
- The reference scales on the instrument handle ensure professional esthetic results.



Delivery forms

OptraSculpt Pad Assortment

- 1 OptraSculpt Pad Instrument
- 1 OptraSculpt Pad Refill/4 mm/60
- 1 OptraSculpt Pad Refill/6 mm/60

OptraSculpt Pad Refill/4 mm/60

OptraSculpt Pad Refill/6 mm/60

OptraSculpt Pad Instrument Refill



This product forms part of our Direct Restoratives category. All the products of this category are optimally coordinated with each other.

Ivoclar Vivadent AG
Bendererstr. 2
9494 Schaan
Liechtenstein
Tel. +423 / 235 35 35
Fax +423 / 235 33 60
www.ivoclarvivadent.com

ivoclar
vivadent
passion vision innovation