

Authorized CAD/CAM **labside** units for Ivoclar Vivadent CAD/CAM materials¹


Overview

Manufacturer of CAD-CAM units	System	IPS e.max [®]						IPS Empress [®] CAD	Tetric [®] CAD	Telio [®] CAD			
		IPS e.max [®] CAD			IPS e.max [®] ZirCAD					Inlays, onlays, veneers and crowns	Inlays, onlays, veneers and crowns	Temporary crowns and bridges ³	Implant-supported hybrid-abutment crowns
		Inlays, onlays, veneers and crowns	Three-unit bridges ²	Implant supported hybrid abutments	Implant-supported hybrid-abutment crowns	Veneering solutions (crowns and bridges)	Crowns						
ivoclar digital [®]	PrograMill [®] One												
	PrograMill [®] PM3/5												
	PrograMill [®] PM7												
ivoclar vivadent [®] passion vision innovation	Zenotec [®] mini												
	Zenotec [®] select												
	Zenotec [®] select hybrid												
AMANNGIRRBACH	ceramill [®] mikro ic												
	ceramill [®] motion 2												
Dentsply Sirona	inLab [®] MC XL												
	inLab [®] MC X5												

¹ The range of products may vary from country to country.





² With up to two connected pontics

³ Up to the second premolar as the terminal abutment



Authorized CAD/CAM **labside** units for Ivoclar Vivadent CAD/CAM materials¹

Overview





		ProArt CAD Wax yellow	ProArt CAD Wax blue	ProArt CAD Splint	ProArt CAD Try-In	ProArt CAD Transfer	ProArt CAD Model
Manufacturer of CAD-CAM units	System	Burnout-restaurations for the press technique	Burnout-restaurations for the casting technique	Occlusal splints, drilling templates	Functional try-in, impression trays / functional trays, bite templates	Transfer templates	Models
	PrograMill® PM3/PM5	●	●	●	●	●	●
	PrograMill® PM7	●	●	●	●	●	●
	Zenotec® mini	●	●	●			●
	Zenotec® select	●	●	●	●	●	●
	Zenotec® select hybrid	●	●	●	●	●	●
	ceramill® motion 2	●					
	ceramill® matik	●					
	inLab® MC X5	●	●	●	●		

¹ The range of products may vary from country to country.

● 98.5 mm discs

Authorized CAD/CAM **labside** units for Ivoclar Vivadent CAD/CAM materials¹

Overview

Manufacturer of CAD-CAM units	System	Ivotion	Ivotion Dent/Multi	Ivotion Base	Colado® CAD CoCr4		Colado® CAD Ti5	
		Full denture	Tooth segments/denture teeth	Denture base	Crowns and bridges (full-contour and frameworks)	Implant supported superstructures	Crowns and bridges (Frameworks)	Implant supported superstructures
	PrograMill® PM3/PM5	●	●	●				
	PrograMill® PM7	●	●	●	●	●	●	●
	Zenotec® mini				●	●		
	Zenotec® select	●*	●	●	●	●		
	Zenotec® select hybrid	●*	●	●	●	●		
	ceramill® motion 2							
	ceramill® matik							
	inLab® MC X5		●	●				

¹ The range of products may vary from country to country.
 * Only in machines with material changer

